

Report R26850

Sample S26289

Andy Colton TESTING
 Test Customer
 Hameldown Road
 Okehampton Business Park
 Okehampton
 Devon EX20 1UB

Purchase Order 111

Telephone 00000000000

Fax

Mobile

Email andy.colton@microanalysis.co.uk


ccEmail lab.elemental@gmail.com

Customer Expected Values

C	60.5
H	12
N	4.5
O	15.0
S	8

Microanalysis Results %	U +/- %
C 60.50 60.40 0.290 %	
H 12.00 12.10 0.215 %	
N 4.50 4.45 0.047 %	
O 14.99 15.01 0.180 %	
S 8.00 8.01 0.112 %	

Average Results Chart



CHNS analyses performed using Dumas combustion method, O analyses performed using Unterzaucher pyrolysis method
 Uncertainty values apply to individual results

Sample description and analyst comments, including any deviations from the standard method

WHITE POWDER
 ANALYSIS AS RECEIVED

Received on 16/01/2024, then analysed on 18/03/2024.

The results above apply only to the sample as received and stated in this report. This report may not be reproduced other than in full without prior written approval of the issuing laboratory



Andrew Colton (Analyst), issued on 18/03/2024